



DECLARATION OF CONFORMITY

Manufacturer: Joinstar Biomedical Technology Co.,Ltd.

Address: 10th Floor, Administration Building, NO.519, XingGuo RD., Yuhang Economic and Technological Development Zone, Hangzhou, Zhejiang, China, 311188

EC Representative's Name: Lotus NL B.V.

EC Representative's Address: Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.

Declares, that the product

Product Name and Model:

COVID-19 Antigen Rapid Test (Latex)
1 Test/Kit, 25 Tests/Kit

as described above are in conformity with the requirements as defined in IVDD98/79/EC Annex III.

Additional information:

Conformity assessment route: Directive 98/79/EC, Annex II

Classification: List Others

I, the undersigned, hereby declare that the medical devices specified above conform with the Directive 98/79/EC on in vitro diagnostic medical devices and pertinent essential requirements.

Date Signed:

2020.09.02

Xuyi ZHOU

General Manager

Joinstar Biomedical Technology Co., Ltd.



Joinstar Biomedical Technology Co., Ltd.

TUV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD
ZERTIFIKAT ◆ CERTIFICATE ◆ 證書 ◆ CERTIFICADO ◆ CERTIFIKAT ◆ CERTIFICATE ◆ 證書 ◆ CERTIFICADO ◆ CERTIFIKAT



Product Service

Certificate

No. Q5 087635 0004 Rev. 01

Holder of Certificate: **JOINSTAR BIOMEDICAL TECHNOLOGY CO., LTD.**
10th Floor, Administration Building
No.519 Xingguo Rd.
Yuhang Economic and Technological Development Zone
311188 Hangzhou
PEOPLE'S REPUBLIC OF CHINA

Facility(ies):
JOINSTAR BIOMEDICAL TECHNOLOGY CO., LTD.
10th Floor, Administration Building, No.519 Xingguo Rd., Yuhang Economic and Technological Development Zone, 311188 Hangzhou, PEOPLE'S REPUBLIC OF CHINA

JOINSTAR BIOMEDICAL TECHNOLOGY CO., LTD.
No. 1 Factory Building, No. 519 Xingguo Rd., Yuhang Economic and Technological Development Zone, 311188 Hangzhou, PEOPLE'S REPUBLIC OF CHINA

Certification Mark:



Scope of Certificate: **Design, Development, Production and Distribution of Biochemical Reagent, ELISA Reagent, Clinical Laboratory Instruments and Rapid Diagnostic Reagents**

Applied Standard(s): EN ISO 13485:2016
Medical devices - Quality management systems - Requirements for regulatory purposes (ISO 13485:2016)
DIN EN ISO 13485:2016

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

Report No.: SH2087401

Valid from: 2020-05-27

Valid until: 2023-05-26

Date, 2020-05-07

Christoph Dicks
Head of Certification/Notified Body



JOINSTAR BIOMEDICAL TECHNOLOGY CO., LTD.
Add: NO.519 Xingguo RD. Yuhang Economic and Technological
Development Zone, Hangzhou, Zhejiang, China
Tel: 0571-89028388 Fax: 0571-89028228
Web: <http://www.joinstar.cn> Email: info@joinstar.cn

Clinical Study (internal clinical assessment) Report

COVID-19 Antigen Rapid Test (Latex)

Prepared By: Yanhua ZHANG Date: 02/09/2020

Reviewed By: Qian XU Date: 02/09/2020

Approval By: Zhong WANG Date: 03/09/2020

1 Purpose:

To determine the performance of COVID-19 Antigen Rapid Test compared to a commercially available competitor assay.

2 Operation Information

2.3 Materials:

COVID-19 Antigen Rapid Lot: COV2008003L-T

Competitor assay: Novel Coronavirus (SARS-COV-2) Real Time Multiplex RT-PCR Kit produced by Shanghai ZJ Bio-Tech Co., Ltd., a commercial SARS-COV-2 kit approved by CFDA, is referred herein as the "gold standard" reagent.

3 Method:

180 samples of the novel coronavirus SARS-CoV-2 and 90 samples of healthy people were included in the testing plan. A total of 90 saliva samples, 90 sputum samples and 90 stool samples were tested in the COVID-19 Antigen Rapid assay and the Novel Coronavirus (SARS-COV-2) real-time multiplex RT-PCR kit, and the clinical sensitivity, specificity and overall agreement between the 2 assays calculated.

4. Clinical research results and analysis

4.1 Clinical research results

4.1.1 The Test results were as follows:

Novel Coronavirus	Covid 19 antigen Rapid Test			RT-PCR
	Saliva	Stool	Sputum	throat swab
Healthy 1	-	-	-	-
Healthy 2	-	-	-	-
Healthy 3	-	-	-	-
Healthy 4	-	-	-	-
Healthy 5	-	-	-	-
Healthy 6	-	-	-	-
Healthy 7	-	-	-	-
Healthy 8	-	-	-	-
Healthy 9	-	-	-	-
Healthy 10	-	-	-	-
Healthy 11	-	-	-	-
Healthy 12	-	-	-	-
Healthy 13	-	-	-	-
Healthy 14	-	-	-	-



JOINSTAR BIOMEDICAL TECHNOLOGY CO., LTD.
 Add: NO.519 Xingguo RD. Yuhang Economic and Technological
 Development Zone, Hangzhou, Zhejiang, China
 Tel: 0571-89028388 Fax: 0571-89028228
 Web: <http://www.joinstar.cn> Email: info@joinstar.cn

Healthy 15	-	-	-	-
Healthy 16	-	-	-	-
Healthy 17	-	-	-	-
Healthy 18	-	-	-	-
Healthy 19	-	-	-	-
Healthy 20	-	-	-	-
Healthy 21	-	-	-	-
Healthy 22	-	-	-	-
Healthy 23	-	-	-	-
Healthy 24	-	-	-	-
Healthy 25	-	-	-	-
Healthy 26	-	-	-	-
Healthy 27	-	-	-	-
Healthy 28	-	-	-	-
Healthy 29	-	-	-	-
Healthy 30	-	-	-	-

DAY	Covid 19 antigen Rapid Test			RT-PCR	
	Saliva	Stool	Sputum	throat swab	throat swab (notes)
Day 1 of disease onset	+	+	+	+	
Day 9 of disease onset	+	+	+	+	
Day 3 of disease onset	+	+	+	+	
Day 6 of disease onset	+	+	+	+	
Day 2 of disease onset	+	+	-	+	
Day 5 of disease onset	+	+	+	+	
Day 8 of disease onset	+	+	+	+	
Day 4 of disease onset	+	+	+	+	
Day 7 of disease onset	+	+	+	+	
Day 7 of disease onset	+	+	+	+	
Day 2 of disease onset	+	+	+	+	
Day 2 of disease onset	-	-	+	+	confirm by retest
Day 2 of disease onset	+	+	+	+	
Day 1 of disease onset	+	+	+	+	
Day 1 of disease onset	+	+	+	+	
Day 6 of disease onset	+	+	+	+	



JOINSTAR BIOMEDICAL TECHNOLOGY CO., LTD.

Add: NO.519 Xingguo RD. Yuhang Economic and Technological Development Zone, Hangzhou, Zhejiang, China

Tel: 0571-89028388 Fax: 0571-89028228

Web: <http://www.joinstar.cn> Email: info@joinstar.cn

Day 4 of disease onset	+	+	+	+	
Day 4 of disease onset	+	+	+	+	
Day 3 of disease onset	+	+	+	+	
Day 3 of disease onset	+	+	+	+	
Day 1 of disease onset	+	+	+	+	
Day 2 of disease onset	+	+	+	+	
Day 3 of disease onset	-	+	+	+	
Day 2 of disease onset	+	+	+	+	
Day 3 of disease onset	+	+	+	+	
Day 6 of disease onset	+	+	+	+	
Day 1 of disease onset	-	-	+	+	confirm by retest
Day 6 of disease onset	+	+	+	+	
Day 8 of disease onset	+	+	+	+	
Day 4 of disease onset	+	+	+	+	
Day 3 of disease onset	+	+	+	+	
Day 2 of disease onset	+	+	+	+	
Day 5 of disease onset	+	+	+	+	
Day 7 of disease onset	+	+	+	+	
Day 9 of disease onset	+	+	+	+	
Day 3 of disease onset	+	+	+	+	
Day 1 of disease onset	-	-	+	+	confirm by retest
Day 2 of disease onset	+	+	+	+	
Day 5 of disease onset	+	+	+	+	
Day 4 of disease onset	+	+	-	+	
Day 3 of disease onset	+	+	+	+	
Day 3 of disease onset	+	+	+	+	
Day 5 of disease onset	+	+	+	+	
Day 5 of disease onset	+	+	+	+	
Day 9 of disease onset	+	+	+	+	
Day 7 of disease onset	+	+	+	+	
Day 2 of disease onset	-	+	+	+	
Day 3 of disease onset	+	+	+	+	
Day 1 of disease onset	-	+	-	+	
Day 1 of disease onset	+	+	+	+	
Day 2 of disease onset	+	+	+	+	
Day 4 of disease onset	+	+	+	+	
Day 3 of disease onset	+	+	+	+	
Day 7 of disease onset	+	+	+	+	

Day 5 of disease onset	+	+	+	+	
Day 6 of disease onset	+	+	+	+	
Day 4 of disease onset	+	+	+	+	
Day 6 of disease onset	+	+	+	+	
Day 1 of disease onset	+	+	+	+	
Day 2 of disease onset	+	+	+	+	

4.2 The above results are summarized as follows:

Saliva Sample		Gold standard reagent		Total
		Positive	Negative	
Test reagent	Positive	54	0	54
	Negative	6	30	36
Total		60	30	90

Sputum Sample		Gold standard reagent		Total
		Positive	Negative	
Test reagent	Positive	57	0	57
	Negative	3	30	33
Total		60	30	90

Stool Sample		Gold standard reagent		Total
		Positive	Negative	
Test reagent	Positive	57	0	57
	Negative	3	30	33
Total		60	30	90

Result analysis

Saliva samples: The Novel Coronavirus Spike Glycoprotein Detection Kit showed 90% sensitivity and 100% specificity in saliva samples.

$$\text{Clinical sensitivity (\%)} = [54 / (54 + 6)] \times 100\% = 90\%$$

$$\text{Clinical specificity (\%)} = [30 / (0 + 30)] \times 100\% = 100\%$$

$$\text{Total agreement rate (\%)} = [(54 + 30) / (54 + 6 + 0 + 30)] \times 100\% = 93.3\%$$

Sputum samples: The Novel Coronavirus Spike Glycoprotein Detection Kit showed 95%



sensitivity and 100% specificity in sputum samples.

Clinical sensitivity (%) = $[57 / (57 + 3)] \times 100\% = 95\%$

Clinical specificity (%) = $[30 / (0 + 30)] \times 100\% = 100\%$

Total agreement rate (%) = $[(57 + 30) / (57 + 3 + 0 + 30)] \times 100\% = 96.7\%$

Stool samples: The Novel Coronavirus Spike Glycoprotein Detection Kit showed 95% sensitivity and 100% specificity in stool samples.

Clinical sensitivity (%) = $[57 / (57 + 3)] \times 100\% = 95\%$

Clinical specificity (%) = $[30 / (0 + 30)] \times 100\% = 100\%$

Total agreement rate (%) = $[(57 + 30) / (57 + 3 + 0 + 30)] \times 100\% = 96.7\%$

5 Discussion and conclusion

The overall clinical performance of the Joinstar COVID-19 Antigen Rapid Test was comparable with the data obtained in the Novel Coronavirus (SARS-COV-2) Real Time Multiplex RT-PCR assay and supports the use of the Joinstar assay in the detection of SARS-CoV-2.